



ORP Platinum series combination electrode, 5-pin, w/ temperature

Product #: 5193700



Hazardous

OBSOLETE ITEM
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

A good performer in dirtier, more difficult samples where redox measurements are commonly made. This platinum-type electrode with temperature sensor features a reference junction that resists clogging and is easy to clean. Includes two gel cartridges and

Accurate

You don't just need an answer, but the right one, and fast. Optimized stabilization algorithms in Hach meter platforms eliminate operator guesswork and reduce inaccurate measurements. Many electrodes work specifically to deliver accurate results and t

Easy to Use

Now, more than ever, you're looking for new ways to stretch your dollar and resources. Skip studying complicated user manuals and spend your time measuring, thanks to Hach's innovative menu designs and simple keypads for truly intuitive operation.

Rugged and Durable

Rely on Hach's durable measuring systems designed to withstand the rigors of outdoor testing and frequent use. Choose from the industry's toughest probes or non-glass electrodes for rugged, waterproof systems that deliver day after day.

Flexible

Need to measure in both the lab and field? Hach offers efficiency with reliable, high-performance systems that easily move with you from the benchtop directly to the testing site. Measure the parameters you need today with the assurance you can expand them tomorrow.

Specifications

Cable Length:	0.91 m
Diameter:	12.3 mm
Electrode Type:	Laboratory
Filling Solution:	2546902
Length:	241 mm
Material:	Sensor Body: PMP
Probe Junction :	Open
Range:	± 2000 mV
Sensor Type :	Platinum annular ring
Special Features:	Platinum series
Temperature Range:	Continuous use: 0 - 45 °C
Temperature Range 2:	N/A N/A N/A
Warranty:	6 months

Required Accessories

- Potassium Chloride Reference Electrolyte Gel Cartridges (Item 2546902)