



## ORP gelfilled combination electrode

Product #:

5193900

 **Hazardous**

**OBSOLETE ITEM**

This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Robust ORP combination probe (non-refillable) for measurement in the field or lab. Temperature range 0 to 80 °C, features 5-pin connector, 1-meter cable (for all Sension meters).

---

### Specifications

Cable Length:	1 m
Diameter:	12 mm
Electrode Type:	Laboratory
Filling Solution:	Non-refillable gel
Length:	135 mm
Material:	Sensor Body: Epoxy
Probe Junction :	Ceramic pin
Range:	± 2000 mV
Sensor Type :	Platinum annular ring
Temperature Range:	Continuous use: 0 - 45 °C
Temperature Range 2:	N/A N/A N/A
Warranty:	6 months