



## Sension Dissolved oxygen probe, 5-pin connector, 3 m cable

Product #:

5197003

 **Hazardous**

### **OBSOLETE ITEM**

This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Clark-Type Polarographic Dissolved Oxygen probe with 3 meter cable and 5-pin circular connector compatible with Hach Sension meters.

### **Accurate**

You don't just need an answer, but the right one, and fast.

### **Easy to Use**

Now, more than ever, you're looking for new ways to stretch your dollar and resources. Skip studying complicated user manuals and spend your time measuring, thanks to Hach's innovative menu designs and simple keypads for truly intuitive operation.

### **Flexible**

Need to measure in both the lab and field? Hach offers efficiency with reliable, high-performance systems that easily move with you from the benchtop directly to the testing site. Measure the parameters you need today with the assurance you can expand them tomorrow.

### **Rugged and Durable**

Rely on Hach's durable measuring systems designed to withstand the rigors of outdoor testing and frequent use. Choose from the industry's toughest probes or non-glass electrodes for rugged, waterproof systems that deliver day after day.

---

## Specifications

Cable Length:	3 m
Electrode Type:	Laboratory Polarographic Dissolved Oxygen
Filling Solution:	2759126
Length:	149 mm
Material:	Sensor Body: Epoxy
Range:	0 - 20 mg/L
Sensor Type :	Platinum cathode, silver anode with replaceable PTFE membrane module (item no. 5197300)
Temperature Range:	0 - 50 °C
Warranty:	6 months

---

## Required Accessories

- Membranes (Item 5197300)
- Filling Solution, Internal, Polarographic DO, 50 mL (Item 2759126)