



## Sension Platinum series combination pH electrode, BNC and DIN connectors

Product #:

5191022



**Hazardous**

**OBSOLETE ITEM**  
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Hach Platinum Series pH probe is a standard combination pH probe with a replaceable gel reference and built-in temperature sensor. Probe comes with two electrolyte cartridges and has a fixed cable with BNC and DIN connectors compatible Hach EC series meters.

### Accurate

You don't just need an answer, but the right one, and fast.

Optimized stabilization algorithms in Hach meter platforms eliminate operator guesswork and reduce inaccurate measurements. Many electrodes work specifically to deliver accurate results and the quickest response, even in challenging environments.

### Easy to Use

Now, more than ever, you're looking for new ways to stretch your dollar and resources. Skip studying complicated user manuals and spend your time measuring, thanks to Hach's innovative menu designs and simple keypads for truly intuitive operation.

### Rugged and Durable

Rely on Hach's durable measuring systems designed to withstand the rigors of outdoor testing and frequent use. Choose from the industry's toughest probes or non-glass electrodes for rugged, waterproof systems that deliver day after day.

### Flexible

Need to measure in both the lab and field? Hach offers efficiency with reliable, high-performance systems that easily move with you from the benchtop directly to the testing site. Measure the parameters you need today with the assurance you can expand in the future.

---

## Specifications

Cable Length:	0.91 m
Connector (text):	BNC and DIN Connectors
Diameter:	12.7 mm
Electrode Type:	Laboratory
Filling Solution:	2546902
Length:	241 mm
Material:	Sensor Body: Epoxy
Probe Junction :	Open
Range:	0 - 14 pH
Sensor Type :	Glass
Temperature Range:	Continuous use: 0 - 45 °C
Temperature Range 2:	N/A N/A N/A
Warranty:	6 months

---

## Required Accessories

- Potassium Chloride Reference Electrolyte Gel Cartridges (Item 2546902)