



## sensION Temperature Probe, DIN connector

Product #:

5198022



**Hazardous**

### **OBSOLETE ITEM**

This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

sensION Temperature Probe has a DIN connector.

### **Accurate**

You don't just need an answer, but the right one, and fast.

### **Easy to Use**

Now, more than ever, you're looking for new ways to stretch your dollar and resources. Skip studying complicated user manuals and spend your time measuring, thanks to Hach's innovative menu designs and simple keypads for truly intuitive operation.

---

## Specifications

Cable Length:	1 m
Connector (text):	DIN Connector
Diameter:	6 mm
Electrode Type:	Laboratory temperature
Length:	120 mm
Material:	Sensor Body: Stainless Steel
Sensor Type :	Thermistor
Temperature Range:	0 - 80 °C
Temperature Range 2:	N/A N/A N/A
Warranty:	6 months
What's included?:	Includes: sensION 51980 Temperature Probe and User Manual

---

## What's included?

Includes: sensION 51980 Temperature Probe and User Manual