



Sension Fluoride ion selective electrode (ISE), BNC connector

Product #: 5192800

OBSOLETE ITEM
This item is no longer available.



Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Combination fluoride ion selective electrode (ISE) with replaceable gel electrolyte and fixed cable with BNC connector. Probe comes with extra replacement electrolyte gel reference cartridge.

Accurate

You don't just need an answer, but the right one, and fast.

Optimized stabilization algorithms in Hach meter platforms eliminate operator guesswork and reduce inaccurate measurements. Many electrodes work specifically to deliver accurate results and the quickest response, even in challenging environments.

Easy to Use

Now, more than ever, you're looking for new ways to stretch your dollar and resources. Skip studying complicated user manuals and spend your time measuring, thanks to Hach's innovative menu designs and simple keypads for truly intuitive operation.

Rugged and Durable

Rely on Hach's durable measuring systems designed to withstand the rigors of outdoor testing and frequent use. Choose from the industry's toughest probes or non-glass electrodes for rugged, waterproof systems that deliver day after day.

Flexible

Need to measure in both the lab and field? Hach offers efficiency with reliable, high-performance systems that easily move with you from the benchtop directly to the testing site. Measure the parameters you need today with the assurance you can expand in the future.

Specifications

Cable Length:	1 m
Diameter:	12 mm
Electrode Type:	Laboratory combination
Filling Solution:	2546902
ISA Required:	258999
Length:	110 mm
Material:	Sensor Body: PMP
Probe Junction :	Open
Range:	0.01 - saturation mg/L F-
Sensor Type :	Crystal membrane
Special Features:	Platinum Series
Temperature Range:	0 - 80 °C
Temperature Range 2:	N/A N/A N/A
Thermistor :	-
Warranty:	6 months

Required Accessories

- Potassium Chloride Reference Electrolyte Gel Cartridges (Item 2546902)
- Fluoride Ionic Strength Adjustor (ISA) Powder Pillows, pk/100 (Item 258999)