



Sension+ pH gel combination electrode (general use)

Product #: LZW5010T.97.002
ZAR Price: Contact Hach
Available

Low cost gel-filled pH electrode with ceramic diaphragm. Built-in Pt1000 for Automatic temperature compensation. Polycarbonate body. Range 0-14 pH, 0-80 °C. 1 m cable. For sensION+ benchtop instruments.

This low maintenance pH electrode is suitable for measurements in aqueous media in general. A storage protector prevents breakage and avoids over pressurisation of the diaphragm.

Automatic temperature compensation

Built-in Pt1000 temperature sensor.

Inclusive storage protector.

Specifications

Accuracy:	0.002 pH with MM374 Meter
Cable Length:	1 m
Connector:	BNC for pH + Banana for Pt 1000
Dimensions (D x L) :	12 mm x 130 mm
Electrode Type:	Laboratory
Filling Solution:	Non-refillable gel
Length:	130 mm
Material:	Sensor Body: Polycarbonate
Parameter:	pH
Probe Junction :	Ceramic pin
Range:	0 - 14 pH
Sensor Type :	Glass
Temperature Range:	0 - 80 °C
Thermistor :	Pt1000
Warranty:	6 months