



## Sension+ pH combination electrode (dirty samples)

**Product #:** LZW5011T.97.002  
**ZAR Price:** Contact Hach  
**Ships within 3-5 weeks**

pH electrode with open circular diaphragm and solid electrolyte. Automatic temperature compensation. Glass body. Range 2-14 pH, 0-80 °C. 1 m cable. For benchtop meter.

The electrode is suitable for waste water, dirty samples, viscous media, and food applications. The open diaphragm around the membrane facilitates the measuring in some media where other electrodes have certain difficulties.

### Automatic temperature compensation

Built-in Pt1000 temperature sensor.

---

## Specifications

Accuracy:	0.002 pH with MM374 Meter
Cable Length:	1 m
Connector:	BNC for pH + Banana for Pt 1000
Dimensions (D x L) :	12 mm x 130 mm
Electrode Type:	Laboratory
Filling Solution:	Non-refillable solid polymer
Length:	130 mm
Material:	Sensor Body: Glass
Parameter:	pH
Probe Junction :	Open
Range:	2 - 14 pH
Sensor Type :	Glass
Special Features:	Solid gel for high solids content
Temperature Range:	0 - 80 °C
Thermistor :	Pt1000
Warranty:	6 months