



Sension+ Multisensor for pH & ORP/Redox

Product #: LZW5045.97.0002
ZAR Price: Contact Hach
Ships within 3-5 weeks

Multisensor for pH (0-14), ORP (± 2000 mV) & temperature (0-80 °C). Automatic temperature compensation. Cable length 1.5m. Specific for Sension+ MM110 multimeter.

With Polycarbonate body, gel electrolyte and ceramic diaphragm.

Attachable calibration tubes can screw into handle of electrode (only 2 mL of buffer or standard). Storage protector prevents breakage and avoids over pressurization of the diaphragm.

Accurate

Easy to use

Rugged & durable

Specifications

Accuracy:	0.02 pH with MM110 Meter
Cable Length:	1.5 m
Connector:	Sension+ MP8
Dimensions (D x L) :	12 mm x 85 mm
Electrode Type:	Portable combination pH, ORP
Filling Solution:	Non-refillable gel
Length:	85 mm
Material:	Sensor Body: Polycarbonate
Parameter:	pH, ORP
Probe Junction :	Ceramic
Range:	ORP: ± 2000 mV
Range 2:	pH: 0 - 14 pH
Sensor Type :	pH: Glass
	ORP: Platinum
Special Features:	pH, ORP
Temperature Range:	0 - 80 °C
Thermistor :	Pt1000

Type:	Portable
Warranty:	6 months