



Sension+ pH gel electrode with ceramic diaphragm (general use)

Product #: LZW5050T.97.002
ZAR Price: Contact Hach
Available

Low cost gel-filled pH electrode with ceramic diaphragm. Automatic temperature compensation. Polycarbonate body. Range 0-14 pH, 0-80 °C. 1 m cable. For portable meters.

This low maintenance pH electrode is suitable for measurements in aqueous media in general. A storage protector prevents breakage and avoids over pressurisation of the diaphragm.

Ergonomic and robust

Integrated temperature sensor

Low maintenance

Specifications

Accuracy:	0.02 pH with pH1 Meter
Cable Length:	1 m
Connector:	Sension+ MP5
Dimensions (D x L) :	12 mm x 85 mm
Electrode Type:	Portable
Filling Solution:	Non-refillable gel
Length:	85 mm
Material:	Sensor Body: Polycarbonate
Parameter:	pH
Probe Junction :	Ceramic pin
Range:	0 - 14 pH
Sensor Type :	Glass
Temperature Range:	0 - 80 °C
Thermistor :	Pt1000
Type:	Portable
Warranty:	6 months