



Sension+ pH combination electrode with open diaphragm (dirty samples)

Product #: LZW5051T.97.002
ZAR Price: Contact Hach
Ships within 3-5 weeks

pH electrode with open circular diaphragm and solid electrolyte. Automatic temperature compensation. Glass body. Range 2-14 pH, 0-80 °C. 1m cable. For portable meters.

The electrode is suitable for (very dirty) waste water, milk, juices, and food applications. A polypropylene working protector prevents accidental breakage of the electrode.

Automatic temperature compensation
Bulit-in Pt1000 temperature sensor.

Ergonomic and robust

Low maintenance

Specifications

Accuracy:	0.02 pH with pH1 Meter
Cable Length:	1 m
Connector:	Sension+ MP5
Dimensions (D x L) :	12 mm x 85 mm
Electrode Type:	Portable
Filling Solution:	Non-refillable solid polymer
Length:	85 mm
Material:	Sensor Body: Glass
Parameter:	pH
Probe Junction :	Open
Range:	2 - 14 pH
Sensor Type :	Glass
Special Features:	With Working Protector. Solid gel for high solids content.
Temperature Range:	0 - 80 °C
Thermistor :	Pt1000
Type:	Portable
Warranty:	6 months