



Sension+ ISE electrode, Sodium (Na⁺), fixed cable 1m, BNC

Product #: LZW9650.97.0002
ZAR Price: Contact Hach
Ships within 3-5 weeks

Ion-selective indicator electrode for Sodium measurements. Range 0.05-23,000 mg/L Na⁺, 0-60 °C. With 1m cable and BNC connector.

Specifications

Accuracy:	±0.02 mV or 0.05%, whichever is greater (application dependent)
Cable Length:	1 m
Connector:	BNC Connector
Dimensions (D x L) :	12 mm x 120 mm
Electrode Type:	Laboratory Half-Cell ISE
ISA Required:	4451569
Length:	120 mm
Material:	Sensor Body: Glass
Parameter:	Sodium
Range:	0.05 - 23000 mg/L Na ⁺
Sensor Type :	Sodium Selective Glass Membrane
Special Features:	Indicator only; requires a reference electrode.
Temperature Range:	0 - 60 °C
Type:	Lab ISE, reference needed
Warranty:	6 months

Required Accessories

- Sension+ Electrolyte Solution 0.1M for ISE Electrodes, 125 mL (Item LZW9901.99)