



Sension+ ISE electrode, Nitrate (NO_3^-), fixed cable 1m, BNC

Product #: LZW9662.97.0002
ZAR Price: Contact Hach
Call for ship date

Ion-selective indicator electrode for Nitrate measurements. Range 0.4-62,000 mg/L NO_3^- , 5-40 °C. With 1m cable and BNC connector.

Specifications

Cable Length:	1 m
Connector:	BNC Connector
Dimensions (D x L) :	12 mm x 120 mm
Electrode Type:	Laboratory Half-Cell ISE
Filling Solution:	-
ISA Required:	2984799
Length:	120 mm
Material:	Sensor Body: Epoxy
Parameter:	Nitrate
Probe Junction :	-
Range:	0.1 mg/L (7×10^{-6} M) - 14,000 mg/L (1 M) NO_3^- -N
Sensor Type :	Solid-state PVC membrane
Special Features:	Dry storage & fast response time. No replacement membranes.
Temperature Range:	0 - 50 °C
Thermistor :	-
Type:	Lab ISE, reference needed
Warranty:	6 months

Required Accessories

- Nitrate Ionic Strength Adjustor Powder Pillows - Intellical, pk/100 (Item 2984799)
- Sension+ Electrolyte Solution 0.1M for ISE Electrodes, 125 mL (Item LZW9901.99)