



Sension+ ISE combination electrode, fluoride, fixed cable 1 m, BNC

Product #: LZW9655C.97.002
ZAR Price: Contact Hach
Ships within 3-5 weeks

Ion-selective combination electrode for Fluoride measurements. Range 0.02-20,000 mg/L F⁻, 5-50 °C. With 1 m cable and BNC connector.

Requires virtually no maintenance

The solid-state sensor design allows for DRY storage of the ISE without a shelf life or membrane replacement.

Provides fast, stable, and accurate response in a variety of sample types

Specifications

Accuracy:	±0.02 mV or 0.05%, whichever is greater (application dependent)
Cable Length:	1 m
Connector:	BNC Connector
Dimensions (D x L) :	12 mm x 120 mm
Electrode Type:	Laboratory Combination ISE
Filling Solution:	Non-refillable driTEK gel
ISA Required:	258999
Length:	120 mm
Material:	Sensor Body: Epoxy
Parameter:	Fluoride
Probe Junction :	Porous PTFE Annular Ring
Range:	0.01 - 19000 mg/L F ⁻
Sensor Type :	Solid-state crystal membrane
Special Features:	Combination electrode. Maintenance free: No electrolyte or spare membranes needed
Temperature Range:	5 - 50 °C
Type:	Lab ISE combination
Warranty:	6 months

Required Accessories

- Fluoride Ionic Strength Adjustor (ISA) Powder Pillows, pk/100 (Item 258999)