



Sension+ 5130 Portable Polarographic Dissolved Oxygen Probe

Product #: LZW5130.97.0002
ZAR Price: Contact Hach
Ships within 3-5 weeks

DO polarographic probe with 1 meter cable and MP5 connector for use with portable meters.

Polarographic dissolved oxygen (DO) probe with a durable plastic body, replaceable oxygen-permeable membrane, refillable reference electrolyte, and built-in temperature sensor. The probe has a fixed 1 meter cable and MP5 connector dedicated for use with Portable Dissolved Oxygen meters. Probe is shipped with 30 mL bottle of electrolyte solution, two membrane modules and polishing strips.

The electrode is ideal for dissolved oxygen measurement in general aqueous applications – particularly in a field environment.

Note: Sension+ 5130 DO Probe is compatible with versions of the Sension+ DO6 meter produced since January 2013. To ensure compatibility with your Sension+ DO6 meter, please contact Support for assistance.

Specifications

Accuracy:	±0.5 mg/L
Cable Length:	1 m
Connector:	Sension+ MP5
Dimensions (D x L) :	12 mm x 120 mm
Electrode Type:	Portable Polarographic Dissolved Oxygen
Filling Solution:	2759123
Length:	120 mm
Material:	Probe body: ABS & POM (nylon)
Parameter:	Dissolved Oxygen
Range:	0.03 mg/L - saturation
Sensor Type :	Replaceable PTFE membrane; Pt cathode; Ag anode
Temperature Range:	0 - 50 °C
Thermistor :	Integrated 30kohm NTC
Type:	Portable
Warranty:	6 months

Required Accessories

- Membranes (Item 5197300)
- Filling Solution, Internal, Polarographic DO, 50 mL (Item 2759126)