



Sension+ 5131 Laboratory polarographic dissolved oxygen probe

Product #: LZW5131.97.0002
ZAR Price: Contact Hach
Available

5131 polarographic DO probe with 3 meter cable and BNC + banana connectors for use with Sension+ MM378 benchtop meter.

Specifications

Accuracy:	±0.5 mg/L
Cable Length:	1 m
Connector:	BNC Connector
Dimensions (D x L) :	12 mm x 120 mm
Electrode Type:	Polarographic Dissolved Oxygen
Filling Solution:	2759123
Length:	120 mm
Material:	Probe body: ABS & POM (nylon)
Parameter:	Dissolved Oxygen
Range:	0.03 mg/L - saturation
Sensor Type :	Replaceable PTFE membrane; Pt cathode; Ag anode
Temperature Range:	0 - 50 °C
Thermistor :	Integrated 30kohm NTC
Warranty:	6 months

Required Accessories

- Membranes (Item 5197300)
- Filling Solution, Internal, Polarographic DO, 50 mL (Item 2759126)