



## Sension+ Chloride analysis package

**Product #:** LZW9652AP.97.02  
**ZAR Price:** Contact Hach  
**Ships within 6-8 weeks**

The Sension+ Chloride ISE Analysis Package contains the Sension+ 9652C combination chloride ISE (LZW9652C.97.002), 100 pack of chloride ISA powders (2318069), 1 L of 100 mg/L chloride standard solution (2370853), 500 mL of 1000 mg/L chloride standard solution (18349).

Sension+ 9652C is a combination chloride ion selective electrode (ISE) with an epoxy body, non-refillable, double junction reference, and non-replaceable crystalline membrane selective to chloride ions in solution. The 9652C has a fixed 1 meter cable and BNC connector for laboratory use with the Hach Sension+ MM340 GLP pH/ISE Laboratory Meter or the Sension+ MM374 Laboratory Multimeter. The 9652C is ideal for measuring chloride concentrations in wastewater, drinking water and general water quality laboratory applications.

### Requires virtually no maintenance

The solid-state sensor design allows for DRY storage of the ISE without a shelf life or membrane replacement.

### Provides fast, stable, and accurate response in a variety of sample types

---

## Specifications

Cable Length:	1 m
Connector:	BNC Connector
Dimensions (D x L) :	12 mm x 120 mm
Length:	120 mm
Parameter:	Chloride
Temperature Range:	5 - 50 °C
Type:	Lab ISE analysis package
Warranty:	6 months
What's included?:	Kit includes: Sension+ 9652C combination Chloride ISE (no. LZW9652C.97.002), two Hach chloride standard solutions (100 mg/L (1 L, no. 2370853) and 1000 mg/L (500 mL, no. 18349)), Hach chloride ionic strength adjustor (ISA) Powder Pillows (pack of 100, no. 2318069), and ISE Basic User Manual.

---

## What's included?

Kit includes: Sension+ 9652C combination Chloride ISE (no. LZW9652C.97.002), two Hach chloride standard solutions (100 mg/L (1 L, no. 2370853) and 1000 mg/L (500 mL, no. 18349)), Hach chloride ionic strength adjustor (ISA) Powder Pillows (pack of 100, no. 2318069), and ISE Basic User Manual.

---

**Required Accessories**

- Chloride Ionic Strength Adjustor (ISA) Powder Pillows, pk/100 (Item 2318069)