



Sension+ pH combination electrode 5012T with A.T.C. for high performance for general aqueous samples

Product #: LZW5012T.97.002

ZAR Price: Contact Hach

High performance pH combination electrode 5012T with ceramic diaphragm, liquid electrolyte (KCl 3M AgCl saturated) and built-in Pt1000 for Automatic Temperature Compensation. Glass body. Range 0-14 pH, 0-80 °C. 1 m cable. For benchtop meter.

The electrode is suitable for aqueous media in general.

Automatic temperature compensation

Built-in Pt1000 temperature sensor.

Specifications

Accuracy:	0.002 pH with MM374 Meter
Cable Length:	1 m
Connector:	BNC for pH + Banana for Pt 1000
Dimensions (D x L) :	12 mm x 130 mm
Electrode Type:	Laboratory
Filling Solution:	LZW9501.99 KCl 3M + AgCl
Length:	130 mm
Material:	Sensor Body: Glass
Parameter:	pH
Probe Junction :	Ceramic
Range:	0 - 14 pH
Sensor Type :	Glass
Temperature Range:	0 - 80 °C
Thermistor :	Pt1000
Warranty:	6 months