



Sension+ pH combination puncture electrode 5053T for food applications, automatic temperature compensation, for portable meters

Product #: LZW5053T.97.002

ZAR Price: Contact Hach

No availability shown

pH puncture electrode with built-in Automatic temperature compensation Pt 1000 for portable meters, solid gel filled, encapsulated Ag/AgCl reference, 1 ceramic and 1 open diaphragm. Glass electrode, 316 stainless steel body. Range 2-14 pH, 0-60 °C. 1m cable, MP5 connector.

The electrode is suitable for food applications, e.g cheese, meat, fish, fruits, bread doughs, etc.

Easy to clean

316 stainless steel body

Automatic temperature compensation

Built-in Pt1000 temperature sensor.

High stability and long lifetime

Encapsulated Ag/AgCl reference system

Specifications

Accuracy:	0.02 pH with pH1 Meter
Cable Length:	1 m
Connector:	Sension+ MP5
Electrode Type:	Portable
Filling Solution:	Non-refillable gel
Material:	External body: 316 stainless steel, pH sensor: Glass
Parameter:	pH
Probe Junction :	1 ceramic, 1 open diaphragm
Range:	2 - 14 pH
Sensor Type :	Glass
Special Features:	Puncture probe, stainless steel body.
Temperature Range:	0 - 60 °C
Thermistor :	Pt1000
Type:	Portable
Warranty:	6 months