



Sension+ 5292 Conductivity cell with two Platinum plates for general measurements in lab.

Product #: LZW5292.97.0002
ZAR Price: Contact Hach

Conductivity cell with two Platinum plate electrodes and glass body for Sension+ benchtop meters. Cell constant $k=1\text{ cm}^{-1}$. 1 m cable.

Measuring range # 100 mS/cm, -10 - 80 °C.

2 pole Platinum plate cell.

Specifications

Cable Length:	1 m
Connector:	Banana (2x)
Dimensions (D x L) :	12 mm x 130 mm
Electrode Type:	Laboratory
Length:	130 mm
Material:	Sensor Body: Glass, Sensor Cell: Platinum
Parameter:	Conductivity
Range:	# 100 mS/cm
Sensor Type :	2-pole Platinum plate, $k = 1\text{ cm}^{-1}$
Temperature Range:	-10 - 80 °C
Warranty:	12 months