



Sension+ 5208 pH liquid filled combination electrode for Micro Sample volumes, 3 mm diameter, screw cap S7

Product #: LZW5208.97.0002
ZAR Price: Contact Hach

pH liquid electrode (refillable) with ceramic diaphragm. Glass body. Range 0-14 pH, 0 to 80 °C. Screw cap S7. 3 mm diameter, 5 mm minimum immersion depth.

This pH electrode is suitable for Micro Sample volumes.

For Micro Sample volumes

High stability reference electrode due to sleeved cartridge

Encapsulated reference system

Specifications

Accuracy:	0.002 pH with pH31 Meter
Connector:	Screw cap S7
Dimensions (D x L) :	3 mm x 70 mm
Electrode Type:	Laboratory
Filling Solution:	LZW9500.99 KCl 3M
Length:	70 mm
Material:	Sensor Body: Glass
Parameter:	pH
Probe Junction :	Ceramic
Range:	0 - 14 pH
Sensor Type :	Glass
Temperature Range:	0 - 80 °C
Warranty:	12 months

Required Accessories

- AS7 cable / 1M / BNC for instruments w/ BNC connector (Item LZW9055.99)