



DR1900 Portable Spectrophotometer with USB + Power Module

Product #:

9624700

ZAR Price:

Contact Hach

No availability shown

Lab methods in the palm of your hand

The DR1900 is perfect for use in the field because it is the lightest and most compact portable spectrophotometer. It also can perform that largest number of tests of any portable instrument. Underneath the rugged exterior, the DR1900 has the largest number--over 220--of the most commonly tested preprogrammed methods already built in. You can also use the easy-to-use interface to create your own methods. Also, it is flexible because you can use the widest range of vial sizes of any other spectrophotometer.

Portability - Easy to hold, and operate, the compact DR1900 is a valuable tool for your field technicians. This is lightest and smallest spec on the market.

Accuracy - Tests are performed with a wavelength range of 340 to 800 nm, which make this a field instrument you can use to find results typically only seen in laboratory instruments.

Ruggedness - The DR1900 is built for rugged conditions. Built as a hand-held spec the DR1900 has a large, clear screen and a simple user interface that makes testing easier than ever in even the most demanding conditions.

Upgrade from your DR2800 or order a new DR1900 today!

Please note: LCK cuvette test evaluation is possible, but it does not have the barcode identification & 10-fold measurement and outlier elimination options available.

Bundle includes the DR1900 Portable Spectrophotometer and USB with power module (LZV813) allowing data transfer and line power options.

Carry anywhere due to a thoughtful, compact design

Test what you need

Ensure accuracy in the field

Make testing easier with clear, easy-to-use interface

Secure in your environment - water and dust-proof design

Specifications

Battery Life: 15 days (typical) @ 5 readings a day / 5 day week without backlight*

* Backlight usage will decrease battery life.

Battery requirements:	4, AA size alkali cells
Cuvette Compatibility 2:	10 mm square / 1 inch square and 13 mm / 16 mm / 1 inch round
Data Logger:	500 measured values (Result, Date, Time, Sample ID, User ID acc. to GLP)
Detector:	Silicon photodiode
Dimensions (H x W x D):	98 mm x 178 mm x 267 mm
Display:	Graphical display 240 x 160 pixel (LCD, b/w, backlit)
Enclosure rating:	IP67
Environmental Conditions: Relative Humidity:	Max. 80 % relative humidity (non-condensing)
Environmental Conditions: Temperature:	10 - 40 °C
Interface:	USB type Mini IP67 (with USB + Power Module)
Mains Connection 2:	110 - 240 V / 50/60 Hz (* optional Module required. May not be available in all regions)
Mains Connection 3:	4x NiMH rechargeable Battery (* optional Module required. May not be available in all region)
Operating conditions:	10 - 40 °C (50 - 104 °F), max. 80 % relative humidity (non-condensing)
Operating Mode:	Transmittance (%), Absorbance and Concentration
Optical System:	Reference beam, spectral
Photometric Accuracy:	± 0.003 Abs @ 0.0 - 0.5 Abs
Photometric Linearity:	< 0.5 % (0.5 - 2.0 Abs)
Photometric Measuring Range:	0 - 3 Abs (wavelength range 340 - 800 nm)
Power Module?:	Yes
Power supply:	4 x AA size Alkaline 4x NiMH rechargeable Battery Power supply: 110 - 240 V; 50/60 Hz
Reproducibility:	± 0.005 Abs (0 - 1 A)
Sample Cell Compatibility:	13 / 16 mm and 1 inch round 10 x 10 mm 1 inch square adapter for round vials and 10x10mm
Source Lamp:	Xenon Flash
Storage conditions:	-30 - 60 °C (-30 - 140 °F), max. 80 % relative humidity (non-condensing)
Stray Light:	< 0.5 %T at 340 nm with NaNO ₂
Temperature Range:	0 - 50 °C (32 - 122 °F)
User Interface Languages:	English, French, German, Italian, Spanish, Portuguese, Bulgarian, Chinese, Czech, Danish, Dutch, Finnish, Greek, Hungarian, Japanese, Korean, Polish, Romanian, Russian, Slovenian, Swedish, Turkish
User programs:	50
Warranty:	12 months
Wavelength Accuracy:	± 2 nm (range 340 - 800nm)
Wavelength Calibration:	Automatic
Wavelength Range:	340 - 800 nm
Wavelength Selection:	Automatic
Weight:	1.5 kg

What's included?:	Includes: portable spectrometer, printed basic instrument manual, dust cover, 1-in. square matched glass sample cells, AA alkaline batteries 4/pk, a four-type vial adapter set, and a USB + Power module (LZV813.99.00002).
-------------------	--

What's included?

Includes: portable spectrometer, printed basic instrument manual, dust cover, 1-in. square matched glass sample cells, AA alkaline batteries 4/pk, a four-type vial adapter set, and a USB + Power module (LZV813.99.00002).