



POCKET Colorimeter II Colorimeter Test Kit for Bromine analysis

Product #: 5870001

 **Hazardous**

OBSOLETE ITEM
This item is no longer available.

Suggested replacements:
• LPV445.97.01110

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Small in size, big on waterproof analysis

Portable Colorimeter programmed for the determination of Bromine, 0.05-4.5 / 0.2-10.0 mg/L Br₂, DPD Method. In carrying case, complete with reagents, cuvettes and operating instructions.

Simple: all functions are available via only four keys

Power: battery operation for a maximum of 2000 tests

Clear readings: even under difficult conditions, thanks to backlit display and large readout

Reliable results with no mains connection

Rugged construction yet lightweight

Specifications

Compliance Certifications:	CE
Detector:	Silicon detector
Dimensions (H x W x D):	155 mm x 61 mm x 35 mm
Display:	LCD, backlit
Enclosure rating:	IP67 (excludes battery compartment)
Lamp Type:	Light Emitting Diode (LED)
Method:	Colorimeter/DPD
Number of tests:	50 - 100
Parameter:	Bromine - As Br ₂
Range:	0.05 - 4.50 mg/L Br ₂
Range 2:	0.2 - 10.0 Br ₂
Smallest increment:	0.01
	0.1
Wavelength:	Fixed wavelength ± 2 nm varies with model
Weight:	0.23 kg