



Be Right™



## POCKET Colorimeter II Colorimeter Test Kit for Aluminium analysis

**Product #:** 5870025  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

### Small in size, big on waterproof analysis

Portable Colorimeter programmed for the determination of Aluminium, 0.1-0.8 mg/L Al, Aluminon Method. In carrying case, complete with reagents, cuvettes and operating instructions.

**Simple: all functions are available via only four keys**

**Power: battery operation for a maximum of 2000 tests**

**Clear readings: even under difficult conditions, thanks to backlit display and large readout**

**Reliable results with no mains connection**

**Rugged construction yet lightweight**

---

### Specifications

Detector:	Silicon detector
Dimensions (H x W x D):	155 mm x 61 mm x 35 mm
Display:	LCD, backlit
Enclosure rating:	IP67 (excludes battery compartment)
Method:	Colorimeter/Aluminon
Number of tests:	100
Parameter:	Aluminum - As Al <sup>3+</sup>
Range:	0.02 - 0.80 mg/L Al
Sample Cell Path Length:	1 cm 25 mm
Smallest increment:	0.01
Warranty:	1 year
Wavelength:	Fixed wavelength ± 2 nm varies with model
Weight:	0.23 kg