



## POCKET Colorimeter II Colorimeter Test Kit for Chlorine Dioxide analysis

Product #: 5870051



**Hazardous**

**OBsolete ITEM**  
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

Suggested replacements:  
• LPV445.97.51110

If necessary, we will change your selected shipping method to accomodate these items.

### Small in size, big on waterproof analysis

Portable Colorimeter programmed for the determination of Chlorine Dioxide, DPD/Glycine Method. In carrying case, complete with reagents, cuvettes and operating instructions.

**Simple:** all functions are available via only four keys

**Power:** battery operation for a maximum of 2000 tests

**Clear readings:** even under difficult conditions, thanks to backlit display and large readout

**Reliable results with no mains connection**

**Rugged construction yet lightweight**

---

### Specifications

Case Style:	D
Compliance Certifications:	CE
Detector:	Silicon detector
Dimensions (H x W x D):	155 mm x 61 mm x 35 mm
Display:	LCD, backlit
Enclosure rating:	IP67 (excludes battery compartment)
Lamp Type:	Lichtemittierende Diode (LED)
Method:	10126
Number of tests:	100
Parameter:	Chlorine Dioxide
Platform:	Powder Pillows
Range:	0.05 - 5.00 mg/L ClO <sub>2</sub>
Smallest increment:	0.01
Wavelength:	Fixed wavelength ± 2 nm varies with model
Weight:	0.23 kg