



Be Right™



## POCKET Colorimeter II Colorimeter Test Kit for Ozone analysis

**Product #:** 5870004  
**ZAR Price:** Contact Hach  
**No availability shown**

### Small in size, big on waterproof analysis

Portable Colorimeter programmed for the determination of Ozone, Indigo-Trisulphate Method. In carrying case, complete with reagents, cuvettes and operating instructions.

**Simple:** all functions are available via only four keys

**Power:** battery operation for a maximum of 2000 tests

**Clear readings:** even under difficult conditions, thanks to backlit display and large readout

**Reliable results with no mains connection**

**Rugged construction yet lightweight**

---

### Specifications

Detector:	Silicon detector
Dimensions (H x W x D):	155 mm x 61 mm x 35 mm
Display:	LCD, backlit
Enclosure rating:	IP67 (excludes battery compartment)
Method:	Colorimeter/ Indigo Trisulfonate
Number of tests:	25
Parameter:	Ozone - As O <sub>3</sub>
Range:	0.01 - 0.25 mg/L O <sub>3</sub>
Range 2:	0.01 - 0.75 mg/L O <sub>3</sub>
Smallest increment:	0.01
Warranty:	1 year
Wavelength:	Fixed wavelength ± 2 nm varies with model
Weight:	0.23 kg