



POCKET Colorimeter II Colorimeter Test Kit for Ozone analysis

Product #: 5870004



Hazardous

OBsolete ITEM
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

Suggested replacements:
• LPV445.97.04110

If necessary, we will change your selected shipping method to accomodate these items.

Small in size, big on waterproof analysis

Portable Colorimeter programmed for the determination of Ozone, Indigo-Trisulphate Method. In carrying case, complete with reagents, cuvettes and operating instructions.

Simple: all functions are available via only four keys

Power: battery operation for a maximum of 2000 tests

Clear readings: even under difficult conditions, thanks to backlit display and large readout

Reliable results with no mains connection

Rugged construction yet lightweight

Specifications

Detector:	Silicon detector
Dimensions (H x W x D):	155 mm x 61 mm x 35 mm
Display:	LCD, backlit
Enclosure rating:	IP67 (excludes battery compartment)
Method:	Colorimeter/ Indigo Trisulfonate
Number of tests:	25
Parameter:	Ozone - As O ₃
Range:	0.01 - 0.25 mg/L O ₃
Range 2:	0.01 - 0.75 mg/L O ₃
Smallest increment:	0.01
Wavelength:	Fixed wavelength ± 2 nm varies with model
Weight:	0.23 kg