



POCKET Colorimeter II Colorimeter Test Kit for Fluoride analysis

Product #:

5870005

 **Hazardous**

OBsolete ITEM
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Small in size, big on waterproof analysis

Portable Colorimeter programmed for the determination of Fluoride, SPADNS Method. In carrying case, complete with reagents, cuvettes and operating instructions.

Simple: all functions are available via only four keys

Power: battery operation for a maximum of 2000 tests

Clear readings: even under difficult conditions, thanks to backlit display and large readout

Reliable results with no mains connection

Rugged construction yet lightweight

Specifications

Detector:	Silicon detector
Dimensions (H x W x D):	155 mm x 61 mm x 35 mm
Display:	LCD, backlit
Enclosure rating:	IP67 (excludes battery compartment)
Method:	8029
Number of tests:	50
Parameter:	Fluoride
Range:	0.1 - 2.0 mg/L
Smallest increment:	0.1
Warranty:	12 months
Wavelength:	Fixed wavelength \pm 2 nm varies with model
Weight:	0.23 kg