



POCKET Colorimeter II Colorimeter Test Kit for Copper analysis

Product #:

5870019



Hazardous

OBsolete ITEM

This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Small in size, big on waterproof analysis

Portable Colorimeter programmed for the determination of Copper, 0.04-5.0 mg/L Cu, Bicinchoninate Method. In carrying case, complete with reagents, cuvettes and operating instructions.

Simple: all functions are available via only four keys

Power: battery operation for a maximum of 2000 tests

Clear readings: even under difficult conditions, thanks to backlit display and large readout

Reliable results with no mains connection

Rugged construction yet lightweight

Specifications

Case Style:	D
Detector:	Silicon detector
Dimensions (H x W x D):	155 mm x 61 mm x 35 mm
Display:	LCD, backlit
Enclosure rating:	IP67 (excludes battery compartment)
Method:	8506
Number of tests:	100
Parameter:	Copper - As free Cu
Range:	0.04 - 5.00 mg/L Cu
Smallest increment:	0.01
Warranty:	12 months
Wavelength:	Fixed wavelength \pm 2 nm varies with model
Weight:	0.23 kg