



## POCKET Colorimeter II Colorimeter Test Kit for Iron analysis

Product #: 5870022

 **Hazardous**

**OBsolete ITEM**  
This item is no longer available.

Suggested replacements:  
• LPV445.97.16110

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

### Small in size, big on waterproof analysis

Portable Colorimeter programmed for the determination of Iron, FerroVer Method. In carrying case, complete with reagents, cuvettes and operating instructions.

**Simple:** all functions are available via only four keys

**Power:** battery operation for a maximum of 2000 tests

**Clear readings:** even under difficult conditions, thanks to backlit display and large readout

**Reliable results with no mains connection**

**Rugged construction yet lightweight**

---

### Specifications

Case Style:	D
Detector:	Silicon detector
Dimensions (H x W x D):	155 mm x 61 mm x 35 mm
Display:	LCD, backlit
Enclosure rating:	IP67 (excludes battery compartment)
Method:	Colorimeter/  FerroVer® (total Fe)
Number of tests:	100
Parameter:	Iron, medium range - As Fe
Range:	0.02 - 5.00 mg/L Fe
Smallest increment:	0.01
Wavelength:	Fixed wavelength $\pm$ 2 nm varies with model
Weight:	0.23 kg