



POCKET Colorimeter II Colorimeter Test Kit for Ammonia analysis

Product #: 5870040



Items with this mark may be considered hazardous under some shipping conditions.

OBsolete ITEM
This item is no longer available.

Suggested replacements:
• LPV445.97.40110

If necessary, we will change your selected shipping method to accomodate these items.

Small in size, big on waterproof analysis

Portable Colorimeter programmed for the determination of Ammonia, 0.01-0.80 mg/L NH₃-N, Salicylate Method. In carrying case, complete with reagents, cuvettes and operating instructions.

Simple: all functions are available via only four keys

Power: battery operation for a maximum of 2000 tests

Clear readings: even under difficult conditions, thanks to backlit display and large readout

Reliable results with no mains connection

Rugged construction yet lightweight

Specifications

Detector:	Silicon detector
Dimensions (H x W x D):	155 mm x 61 mm x 35 mm
Display:	LCD, backlit
Enclosure rating:	IP67 (excludes battery compartment)
Method:	8155
Number of tests:	100
Parameter:	Nitrogen, Ammonia, mid range - For freshwater or seawater**
Range:	0.01 - 0.80 mg/L NH ₃ -N
Sample Cell Path Length:	1 cm 25 mm
Smallest increment:	0.1
Wavelength:	Fixed wavelength ± 2 nm varies with model
Weight:	0.23 kg