



## POCKET Colorimeter II Colorimeter Test Kit for Sulphate analysis

Product #:

5870029



**Hazardous**

**OBsolete ITEM**

This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

### Small in size, big on waterproof analysis

Portable Colorimeter programmed for the determination of Sulphate, 2-70 mg/L  $\text{SO}_4$ , Bariumsulphate Method. In carrying case, complete with reagents, cuvettes and operating instructions.

**Simple:** all functions are available via only four keys

**Power:** battery operation for a maximum of 2000 tests

**Clear readings:** even under difficult conditions, thanks to backlit display and large readout

**Reliable results with no mains connection**

**Rugged construction yet lightweight**

---

### Specifications

Case Style:	D
Detector:	Silicon detector
Dimensions (H x W x D):	155 mm x 61 mm x 35 mm
Display:	LCD, backlit
Enclosure rating:	IP67 (excludes battery compartment)
Method:	8051
Number of tests:	100
Parameter:	Sulfate - As $\text{SO}_4^{-2}$
Range:	2 - 70 mg/L
Smallest increment:	1
Warranty:	12 months
Wavelength:	Fixed wavelength $\pm 2$ nm varies with model
Weight:	0.23 kg