



POCKET Colorimeter II Colorimeter Test Kit for Dissolved Oxygen (DO)

Product #:

5870003

 **Hazardous**

OBsolete ITEM
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

Suggested replacements:

- LPV445.97.03110

If necessary, we will change your selected shipping method to accomodate these items.

Portable Colorimeter programmed for the determination of Dissolved Oxygen, HRDO method. In carrying case, complete with reagents, cuvettes and operating instructions.

Specifications

Compliance Certifications:	CE
Detector:	Silicon detector
Dimensions (H x W x D):	155 mm x 61 mm x 35 mm
Display:	LCD, backlit
Enclosure rating:	IP67
Lamp Type:	Light Emitting Diode (LED)
Method:	Colorimeter/ AccuVac®
Number of tests:	25
Parameter:	Dissolved Oxygen - As DO
Range:	0.2 - 10.0 mg/L O ₂
Smallest increment:	0.1
Wavelength:	Fixed wavelength ± 2 nm varies with model
Weight:	0.23 kg