



Pocket Colorimeter II, CL/pH, replacement instrument

Product #:

5953012

 **Hazardous**

OBSOLETE ITEM

This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

This product is available for purchase solely as a replacement.

*This product has not been evaluated to test for chlorine and chloramines in medical applications in the United States.

Measures one or two parameters

Lightweight and battery operated, suitable for extended field work or quick, on-the-spot process monitoring. Measurement can be made in two ranges for a number of important parameters.

Power: battery operation for a maximum of 2000 tests

The Pocket Colorimeter II is simple enough to use that anyone can obtain reliable results with accuracy comparable to that of a lab instrument.

Factory-Programmed

Each Pocket Colorimeter II Instrument is pre-programmed for one of more than 35 parameters. Many are based on EPA-approved methods. No manual calibration is ever required. Simply zero the instrument with a blank, insert the reacted sample, and read the result.

Waterproof, Lightweight, Rugged

The Pocket Colorimeter II Instrument really can be carried in your pocket! Weighing in at only .23 kg (8.1 oz), the tough, impact-resistant shell protects electronics and optics for years of trouble-free operation. IP67 rating: the Pocket Colorimeter II can remain submerged at 1 meter for 30 minutes and still function properly. And, the instrument is designed to float!

Backlit Display, Data Logging, Recall

Clear readings even under difficult conditions, thanks to lighting and large figures on the display. Record 10 most recent data points along with time of measurement. No more manual record keeping hassles with carrying paper and pen.

Specifications

Dimensions (H x W x D):	32 mm x 62 mm x 157 mm
Enclosure rating:	IP67 (excludes battery compartment)
Parameter:	Chlorine, pH
Wavelength:	± varies with model
Weight:	0.23 kg
What's included?:	Instrument only

What's included?

