



POCKET Colorimeter II Colorimeter Test Kit for Zinc analysis

Product #: 5870009



OBsolete ITEM
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

Suggested replacements:
• LPV445.97.09110

If necessary, we will change your selected shipping method to accomodate these items.

Portable Colorimeter programmed for the determination of Zinc, Zincon method. In carrying case, complete with reagents, cuvettes and operating instructions.

Simple: all functions are available via only four keys

Power: battery operation for a maximum of 2000 tests

Clear readings: even under difficult conditions, thanks to backlit display and large readout

Reliable results with no mains connection

Rugged construction yet lightweight

Specifications

Case Style:	R
Compliance Certifications:	European CE mark
Detector:	Silicon detector
Dimensions (H x W x D):	155 mm x 61 mm x 35 mm
Display:	LCD, backlit
Enclosure rating:	IP67
Lamp Type:	Light Emitting Diode (LED)
Method:	Colorimeter/Zincon
Number of tests:	100
Parameter:	Zinc - As Zn
Range:	0.02 - 3.00 mg/L Zn
Smallest increment:	0.01
Wavelength:	Fixed wavelength \pm 2 nm varies with model
Weight:	0.23 kg