



POCKET Colorimeter II Colorimeter Test Kit for Nickel and Cobalt analysis

Product #:

5870020

 **Hazardous**

OBSOLETE ITEM
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Portable Colorimeter programmed for the determination of Nickel & Cobalt, PAN method. In carrying case, complete with reagents, cuvettes and operating instructions.

Simple: all functions are available via only four keys

Power: battery operation for a maximum of 2000 tests

Clear readings: even under difficult conditions, thanks to backlit display and large readout

Reliable results with no mains connection

Rugged construction yet lightweight

Specifications

Case Style:	D
Compliance Certifications:	European CE mark
Detector:	Silicon detector
Dimensions (H x W x D):	155 mm x 61 mm x 35 mm
Display:	LCD, backlit
Enclosure rating:	IP67
Lamp Type:	Light Emitting Diode (LED)
Method:	8150
Number of tests:	100
Parameter:	Cobalt and Nickel
Range:	0.01 - 1.00 mg/L Ni / 0.02 - 2.00 mg/L Co
Range 2:	0.02 - 2.00 mg/L Co
Sample Cell Path Length:	1 cm 25 mm
Smallest increment:	0.01
	0.01
Warranty:	12 months
Wavelength:	Fixed wavelength \pm 2 nm varies with model

Weight:

0.23 kg