



Test kit TPH in soil, pocket colorimeter II

Product #:

2775000

 **Hazardous**

OBSOLETE ITEM

This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Specifications

Case Style:	Q (2 cases)
Compliance Certifications:	European CE mark
Detector:	Silicon detector
Dimensions (H x W x D):	155 mm x 61 mm x 35 mm
Display:	LCD, backlit
Enclosure rating:	IP67
Instrument:	Kit includes PC II
Lamp Type:	Light Emitting Diode (LED)
Method:	10050
Number of tests:	Up to 18
Parameter:	TPH in Soil**
Range:	20, 50, 100, 200 ppm thresholds
Test Kit Model:	Pocket Colorimeter™ II**
Warranty:	12 months
Wavelength Accuracy:	Fixed wavelength ± 2 nm varies with model
Weight:	0.23 kg