



TPH in Water, Test Kit

Product #: 2774200



Hazardous

OBSOLETE ITEM
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For personnel responsible for underground storage tanks, off-site waste impounds, petrochemical plants, and waste remediation sites, Hach TPH kits provide fast accurate results at low levels. This water kit has thresholds of 2, 5, 10, and 20 ppm.

Provides fast, accurate results at low levels

Pocket Colorimeter II is waterproof, IP67 rated, and datalogs last 10 points with measurement time

Used by personnel responsible for underground storage tanks, off-site waste impounds, petrochemical plants, and waste remediate sites

Specifications

| | |
|----------------------------|---|
| Case Style: | Q |
| Compliance Certifications: | European CE mark |
| Detector: | Silicon detector |
| Dimensions (H x W x D): | 155 mm x 61 mm x 35 mm |
| Display: | LCD, backlit |
| Enclosure rating: | IP67 |
| Instrument: | Kit includes PCII |
| Lamp Type: | Light Emitting Diode (LED) |
| Method: | 10050 |
| Number of tests: | Up to 18 |
| Parameter: | TPH in Water |
| Range: | 2, 5, 10, 20 ppm |
| | thresholds |
| Test Kit Model: | Pocket |
| | Colorimeter™ II |
| Warranty: | 12 months |
| Wavelength Accuracy: | Fixed wavelength \pm 2 nm varies with model |
| Weight: | 0.23 kg |

| | |
|-------------------|--|
| What's included?: | Includes Pocket Colorimeter II, TenSette Pipet, , reagents, Laboratory Marker, Talking Timer, Disposable Wipes, Deionized Water, instruction sheet, and carrying case. |
|-------------------|--|

What's included?

Includes Pocket Colorimeter II, TenSette Pipet, , reagents, Laboratory Marker, Talking Timer, Disposable Wipes, Deionized Water, instruction sheet, and carrying case.