



Conductivity Standard Solution, 18000 µS/cm, NaCl, 100 mL

Product #: 2307442
ZAR Price: Contact Hach
Ships within 3 days

Conductivity standard solutions for laboratory and field analyses

Conductivity standard solution, 18000 µS/cm \pm 50 µS/cm at 25°C, 10246 mg/L \pm 100 mg/L as Sodium Chloride (NaCl), (1 µS/cm \sim 0.5 mg/L TDS), 100 mL bottle.

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate conductivity measurements is good care and calibration of the conductivity cell.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifications

Concentration:	18000 µS/cm at 25°C \pm 50 µS/cm at 25°C
Concentration 2:	10246 mg/L \pm 100 mg/L as Sodium Chloride (NaCl)
Description:	Conductivity calibration solution
Formula:	NaCl
NIST Traceable:	Yes
Package Type:	Bottle
Parameter:	Conductivity
Platform:	Solution
Quantity:	100 mL
Special Remarks:	Observed shelf life after opening, following recommended storage and use: 3 months
Storage conditions:	10 - 25 °C (protect from light)