



## Sodium chloride standard solution, 100 $\mu\text{S}/\text{cm}$ (NIST), 500 mL

**Product #:** 2971849  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

Conductivity standard solution, 100  $\mu\text{S}/\text{cm}$   $\pm 1$   $\mu\text{S}/\text{cm}$  at 25°C, 49 mg/L  $\pm 0.25$  mg/L as Sodium Chloride (NaCl), (1  $\mu\text{S}/\text{cm}$   $\sim$  0.5 mg/L TDS), 500 mL bottle.

---

### Specifications

Concentration:	100 $\mu\text{S}/\text{cm}$ at 25°C $\pm 1$ $\mu\text{S}/\text{cm}$ at 25°C
Concentration 2:	49 mg/L $\pm 0.25$ mg/L as Sodium Chloride (NaCl)
Description:	Conductivity calibration solution
Package Type:	Bottle
Parameter:	Conductivity
Platform:	Solution
Quantity:	500 mL
Special Remarks:	Observed shelf life after opening, following recommended storage and use: 3 months
Storage conditions:	10 - 25 °C (protect from light)