



Natural Water Standard Solution, 3000 ppm Total Dissolved Solids (TDS), 50 mL

Product #: 2974826
ZAR Price: Contact Hach
Ships within 3 days

Conductivity standard solutions for laboratory and field analyses

Natural Water Standard Solution. 40% sodium sulfate/40% sodium bicarbonate/20% sodium chloride. TDS: 3000 ppm \pm 30 ppm, Conductivity: 3960 μ S/cm \pm 39.6 μ S/cm at 25°C. 50 mL pour-top bottle.

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate conductivity measurements is good care and calibration of the conductivity cell.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifications

Concentration:	3000 ppm \pm 30 ppm
Concentration 2:	3960 μ S/cm at 25°C \pm 39.6 μ S/cm at 25°C
Description:	TDS/Conductivity calibration solution
Formula:	Na ₂ SO ₄ , Na ₂ CO ₃ , NaCl
NIST Traceable:	Yes
Package Type:	Bottle
Parameter:	TDS
Platform:	Solution
Quantity:	50 mL
Special Remarks:	Observed shelf life after opening, following recommended storage and use: 3 months
Storage conditions:	10 - 25 °C (protect from light)