



Conductivity Standard Solution, 10000 $\mu\text{S}/\text{cm}$, NaCl, 500 mL

Product #: 2972249
ZAR Price: Contact Hach
Ships within 3 days

Conductivity standard solutions for laboratory and field analyses

Conductivity standard solution, 10000 $\mu\text{S}/\text{cm}$ ± 100 $\mu\text{S}/\text{cm}$ at 25°C, 5.556 g/L ± 0.03 g/L as Sodium Chloride (NaCl), (1 $\mu\text{S}/\text{cm}$ \sim 0.5 mg/L TDS), 500 mL bottle.

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate conductivity measurements is good care and calibration of the conductivity cell.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifications

Concentration:	10000 $\mu\text{S}/\text{cm}$ at 25°C ± 100 $\mu\text{S}/\text{cm}$ at 25°C
Concentration 2:	5.556 g/L ± 0.03 g/L NaCl
Description:	Conductivity calibration solution
Package Type:	Bottle
Parameter:	Conductivity
Platform:	Solution
Quantity:	500 mL
Storage conditions:	10 - 25 °C (protect from light)