



Singlet Single-use Conductivity Standard Solution, 1413 $\mu\text{S}/\text{cm}$, KCl, 20 mL

Product #: 2771420
ZAR Price: Contact Hach
Ships within 3-5 weeks

Singlet Solution Packs

Conductivity Standard Solution. Sealed foil pouches of 1413 $\mu\text{S}/\text{cm} \pm 0.1\%$ at 25 °C. Potassium Chloride (KCl). Ready for use, 20 mL each

Ideal for on-site analyses

Sealed pouches are easy to transport - no risk of spillage.

Easy to Use

Just open the pouch and insert the electrode.

Fresh and reliable

Foil pouches prevent contamination and air exposure.

Design for use with Hach meters portfolio and other meters.

The basis for accurate Conductivity measurements is good care and calibration of the Conductivity cell.

Specifications

Concentration:	1413 $\mu\text{S}/\text{cm}$ at 25°C $\pm 0.1\%$ at 25°C
Concentration 2:	NA
Description:	Conductivity calibration solution
Formula:	KCl
NIST Traceable:	Yes
Package Type:	Singlet
Parameter:	Conductivity
Quantity:	20 mL
Special Remarks:	Discard solution after use
Storage conditions:	10 - 25 °C (protect from light)