



## Singlet Single-use Conductivity Standard Solution, 12.88 mS/cm, KCl, 20 mL

**Product #:** 2771520  
**ZAR Price:** Contact Hach  
**Ships within 3-5 weeks**

### Singlet Solution Packs

Conductivity Standard Solution. Sealed foil pouches of 12.88 mS/cm  $\pm 10\%$  at 25 °C. Potassium Chloride (KCl). Ready for use, 20 mL each

#### Ideal for on-site analyses

Sealed pouches are easy to transport - no risk of spillage.

#### Easy to Use

Just open the pouch and insert the electrode.

#### Fresh and reliable

Foil pouches prevent contamination and air exposure.

#### Design for use with Hach meters portfolio and other meters.

The basis for accurate Conductivity measurements is good care and calibration of the Conductivity cell.

---

### Specifications

Concentration:	12.88 mS/cm at 25°C $\pm 10\%$ at 25°C
Concentration 2:	NA
Description:	Conductivity calibration solution
Formula:	KCl
NIST Traceable:	Yes
Package Type:	Singlet
Parameter:	Conductivity
Quantity:	20 mL
Special Remarks:	Discard solution after use
Storage conditions:	10 - 25 °C (protect from light)