



Ammonia Ionic Strength Adjustor (ISA) Powder Pillows, pk/100

Product #:

4447169

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Ionic Strength Adjustor (ISA) for Ammonia determinations by Ion Selective Electrode (ISE) method. Powder form prevent from diluting samples. Ideal for low Ammonia level readings. Eliminate most interferences. Pack of 100 powder pillows.

Specifications

Concentration:	ISA
Package Type:	Powder Pillows
Parameter:	Ammonium
Quantity:	100 /pk