



Nitrate Ionic Strength Adjustor (ISA) Powder Pillows, pk/100

Product #:

4456369

ZAR Price:

Contact Hach

Ships within 3 days

Ionic Strength Adjustor (ISA) for Nitrate determinations by Ion Selective Electrode (ISE) method. Powder form prevent from diluting samples. Ideal for low Nitrate level readings. Eliminate most interferences. Pack of 100 powder pillows.

Specifications

Concentration:	ISA
Description:	Nitrate Ionic Strength Adjuster Powder Pillows
Method:	Ion Selective Electrode
Number of tests:	100
Package Type:	Powder Pillows
Parameter:	Nitrate
Quantity:	100 /pk
Range:	0.05 - 14000 mg/L NO ₃ ⁻ -N