



Electrode Cleaning Solution for Regular Maintenance, 500 mL

Product #:

2965249

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Electrode cleaning solution for pH electrode regular cleaning maintenance. Convenient for soft cleaning of the pH sensing element after use in clean water samples. 500 mL bottle.

Specifications

Description:	Regular cleaning maintenance
Formula:	KATHON CG, DECONEX 11
Platform:	Solution
Quantity:	500 mL
Solution Type:	Cleaning