



## Nitrate Ionic Strength Adjustor (ISA) Buffer Solution, 500 mL

|                   |                |
|-------------------|----------------|
| <b>Product #:</b> | <b>2488349</b> |
| ZAR Price:        | Contact Hach   |
| Available         |                |

For sample treatment in nitrate determination by ion selective electrode. 500 mL bottle.

---

### Specifications

|              |  |
|--------------|--|
| Description: | Nitrate Ionic Strength Adjuster Solution, 500 mL |
| Method:      | Ion Selective Electrode                          |
| Platform:    | Solution   |