



Ammonia Ionic Strength Adjustor (ISA), Solution, 500 mL

Product #:

2824349

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Ionic Strength Adjustor (ISA) for Ammonia determinations by Ion Selective Electrode (ISE) method. Eliminate most interferences. Solution, 500 mL bottle.

Specifications

Concentration:	ISA
Package Type:	Bottle
Parameter:	Ammonium
Quantity:	500 mL