



## Buffer Powder Pillows, pH 7.00 Colour-coded Yellow, pk/50

**Product #:** 2227066  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

### pH Buffer standard solutions for laboratory and field analyses

pH Buffer Powder Pillows for pH electrode calibration. Each powder pillow prepares 50 mL of fresh buffer solution. Colour-coded for easy identification (Yellow). pH is  $7.00 \pm 0.02$  at 25°C. Pack of 50 Powder Pillows

#### Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard pH Buffer Solutions that are used are key.

#### Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

#### Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

---

### Specifications

Color coded?:	Yes
Description:	pH Buffer Calibration reagent
NIST Traceable:	Yes
Package Type:	Powder Pillows
pH Value(s):	$7.00 \pm 0.02$ pH at 25°C
Quantity:	50 pcs
Special Remarks:	Discard prepared solution after use
Storage conditions:	10 - 25 °C (protect from light)

---

### Required Accessories

- Beaker for pH Electrode Calibration, Yellow, 30 mL, pk/80 (Item SM5012)
- Clippers (Shears) (Item 2369400)
- Wipes, Disposable, 11 x 22 cm, 280/pk (Item 2096900)
- Deionized water, 500 mL (Item 27248)